




# EXTERNAL FIXATION TGF

## ВНЕШНИЙ ФИКСАТОР TGF

Назад, к основам простоты,  
стабильности и универсальности


*and versatility*



IT'S GOT TO MAKE SENSE





## ГЕХФИКС СИСТЕМА НАТЯЖНОЙ ФИКСАЦИИ ПЕРЕЛОМОВ ТГФ



### ФИКСАТОР ТГФ: LE FIXATOR :



- набор для лечения переломов верхних конечностей, состоящий из двух проводников спиц для натяжения, необходимого при репозиции
- с момента начала использования системы ТГФ мы больше не нуждаемся в других методах лечения хирургических переломов шейки плечевой кости
- при помощи этого метода и этой системы стало возможным лечение переломов шейки плечевой кости простым и доступным способом



1 Попытка репозиции посредством внешней манипуляции

2 Использование 5-миллиметровой спицы для репозиции и введение проводника

3 Установка первого проводника спиц Киршнера




4 Установка второго проводника спиц Киршнера для управления поворотом

5 Установка самосверлящей спицы после просверливания диафиза

6 Натяжение проводников спиц Киршнера и установка системы ТГФ





**Connection**

diameter 9mm/9mm

1010-1040

**S**

**Reduction Sleeve**

diameter 9mm/2,5mm

1010-1055

diameter 9mm/4.0mm

1010-1053

diameter 9mm/5mm

1010-1054

**L**

**Double-L Bar**

left rod diameter 5mm

1010-1081

right rod diameter 5mm

1010-1082



**Carbon rod**

diameter 9mm, length 110mm

1010-1064

diameter 9mm, length 150mm

1010-1065



**Lengthening /compression tube**

1010-1100



**Kirschner wire**

diameter 2,5mm, length 300mm

1020-1032

diameter 2,5mm, length 300mm, 25 mm thread middle

1020-1033

diameter 2,5mm, length 250mm tip threaded

1020-1034



**Self drilling /self tapping pin**

Pin diameter 4mm, length 100mm, thread length 40mm

2015-1020

Pin diameter 4mm, length 145mm, thread length 45mm

2015-1025

Pin diameter 4mm, length 180mm, thread length 50mm

2015-1030

Pin diameter 5mm, length 210mm, thread length 55mm

2020-1030



**Trocar and sleeve**

diameter 5mm, length 70mm

5025-1000

diameter 5mm, length 120mm

5025-1200



**Handle for trocar and sleeve**

5025-1100



**Screw handle**

5015-1000



**Drill brace**

5010-1000



**Sterilisation Tray**

6010-1007

Distributed by :



Gexfix S.à.  
Av. De La Praille 50  
CH-1227 Carouge  
Tel: +41 22 308 16 70  
Fax: +41 22 308 16 79  
www.gexfix.ch  
ISO 9001/2000  
ISO 13485/2003

© Eo4sə  
EN 93/42  
FDA



Cat. No: TGF 05-2006

To ensure the best quality of its products and their improvements,  
Gexfix reserves the right to modify all or part of their products